



eva

CASE HISTORIES

VULVAR REJUVENATION

CONDITION: loose vulvar tissue



BEFORE



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 8-12%



VULVAR REJUVENATION & VVA

CONDITION: loose vulvar tissue



BEFORE



DURING



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation + Vulvo Vaginal Atrophy

SESSION DURATION: VER 10 min. 15 min.

POWER: VER 8-12% | VVA 20-22%

SESSION 3

TREATMENT: Vulvar External Rejuvenation + Vulvo Vaginal Atrophy

SESSION DURATION: VER 10 min. | VVA 15 min.

POWER: VER 10-12% | VVA 24-26%

VULVAR REJUVENATION

CONDITION: loose vulvar tissue



BEFORE



AFTER

SESSION 4

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 8-13%



VULVAR REJUVENATION

CONDITION: loose vulvar tissue



BEFORE




AFTER

SESSION 3

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 9-14%



VULVAR REJUVENATION

CONDITION: loose vulvar tissue



BEFORE



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 11%

VULVO VAGINAL ATROPHY

CONDITION: vulvar and vaginal atrophy due to pharmacology



BEFORE



AFTER

SESSION 1

TREATMENT: Vulvo-Vaginal Atrophy

SESSION DURATION: 10 min.

POWER: 15 - 20%



VULVAR REJUVENATION

CONDITION: vulvar lichen



BEFORE



DURING



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation - vulvar tip

SESSION DURATION: 10 min.

POWER: 10%

VULVAR REJUVENATION

CONDITION: vulvar vestibulitis



BEFORE



AFTER

SESSION 3

TREATMENT: Vulvar External Rejuvenation - vulvar tip

SESSION DURATION: 10 min.

POWER: 10%



VULVAR REJUVENATION

CONDITION: vulvar lichen



BEFORE



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation + UPR (hyaluronic acid) - vulvar tip

SESSION DURATION: 10 min.

POWER: 10%

VULVAR REJUVENATION

CONDITION: lichen vulvar



BEFORE



AFTER

SESSION 3

TREATMENT: Vulvar External Rejuvenation + UPR (hyaluronic acid) - vulvar tip

SESSION DURATION: 10 min.

POWER: 10%

VAGINAL LAXITY & VULVAR REJUVENATION

CONDITION: post partum bladder prolaps



BEFORE



AFTER

SESSION 4

TREATMENT: Vaginal Laxity + UPR (hyaluronic acid) - vaginal tip
+ pelvic floor exercise at home

SESSION DURATION: 12 min.

POWER: 20-25%

VULVAR REJUVENATION

CONDITION: vulvar lichen



BEFORE



AFTER

SESSION 3

TREATMENT: Vulvar External Rejuvenation + UPR

SESSION DURATION: 10 min.

POWER: 10%

VULVAR REJUVENATION

CONDITION: lichen vulvar



BEFORE



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 10%

VULVAR REJUVENATION

CONDITION: lichen in the young patient



BEFORE



AFTER

SESSION 2

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 10%



VULVAR REJUVENATION

CONDITION: vulvar atrophy



BEFORE



AFTER

SESSION 3

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 10-12%



VULVAR REJUVENATION & VVA

CONDITION: vulvo vaginal atrophy



BEFORE



AFTER

SESSION 2

TREATMENT: Vulvar External Rejuvenation + Vulvo Vaginal Atrophy

SESSION DURATION: 10 min.

POWER: VER 10-12% VVA 24-26%

VULVAR REJUVENATION

CONDITION: vulvo vaginal atrophy



BEFORE



AFTER

SESSION 5

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: VER 12-15% VVA 26%

VULVAR REJUVENATION

CONDITION: hypotrophy vulvar



BEFORE



AFTER

SESSION 3

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 10%



VULVAR REJUVENATION

CONDITION: vulvar varices



BEFORE




AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation

SESSION DURATION: 10 min.

POWER: 10%



VULVAR REJUVENATION & LASSITY

CONDITION: labia minora hypertrophy



BEFORE



AFTER

SESSION 2

TREATMENT: Vulvar External Rejuvenation + Vaginal Lassity

SESSION DURATION: 10 min.

POWER: VER 10%

VL 25-29%



VULVAR REJUVENATION

CONDITION: laxity



BEFORE



AFTER

SESSION 1

TREATMENT: Vulvar External Rejuvenation + Vaginal Laxity

SESSION DURATION: 10 min.

POWER: VER 10%

VL 25-29%